

ABX Minocal



LOT

CX 463

CAL



(Exp.) 2021-11-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40	4.40							± 0.06
HB HGB	g/dl	12.8	12.8							± 0.2
	g/l	128	128							± 2
	mmol/l	7.95	7.95							± 0.12
HT HCT	%	37.0	37.0							± 1.0
	l/l	0.370	0.370							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	250	250							± 10
VPM MPV	µm ³ ; fl	9.4	9.4							± 0.5

Ref: TEMP-0387 Rev.47 BACK / VERSO 1300060881