

ABX Minocal



LOT

CX 433

CAL



(Exp.) 2019-05-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.85	8.85							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.42							± 0.06
HB HGB	g/dl	12.9	12.9							± 0.2
	g/l	129	129							± 2
	mmol/l	8.01	8.01							± 0.12
HT HCT	%	37.1	36.9							± 1.0
	l/l	0.371	0.369							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	259	255							± 10
VPM MPV	µm ³ ; fl	11.1	10.1							± 0.5

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060881