

# ABX Minocal



**LOT**

**CX 403**

**CAL**



**(Exp.) 2016-11-05**

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H500 OT v1 only	H500 OT H500 CT							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.7</b>	<b>8.8</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.53</b>	<b>4.55</b>							± 0.06
HB HGB	g/dl	<b>13.2</b>	<b>13.3</b>							± 0.2
	g/l	<b>132</b>	<b>133</b>							± 2
	mmol/l	<b>8.20</b>	<b>8.26</b>							± 0.12
HT HCT	%	<b>37.7</b>	<b>38.4</b>							± 1.0
	l/l	<b>0.377</b>	<b>0.384</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>252</b>	<b>260</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>10.0</b>	<b>9.6</b>							± 0.5

Ref: TEMP-0387 Rev.42 BACK / VERSO 1300019741-A