

ABX Minocal



LOT

CX 403

CAL



(Exp.) 2016-11-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H500 OT v1 only	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.7	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.55							± 0.06
HB HGB	g/dl	13.2	13.3							± 0.2
	g/l	132	133							± 2
	mmol/l	8.20	8.26							± 0.12
HT HCT	%	37.7	38.4							± 1.0
	l/l	0.377	0.384							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	252	260							± 10
VPM MPV	µm ³ ; fl	10.0	9.6							± 0.5

Ref: TEMP-0387 Rev.42 BACK / VERSO 1300019741-A