

# ABX Minocal



**LOT**

**CX 460**

**CAL**



**(Exp.) 2021-08-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.70</b>	<b>8.70</b>							± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.55</b>	<b>4.55</b>							± 0.06
HB HGB	g/dl	<b>13.4</b>	<b>13.4</b>							± 0.2
	g/l	<b>134</b>	<b>134</b>							± 2
	mmol/l	<b>8.32</b>	<b>8.32</b>							± 0.12
HT HCT	%	<b>38.7</b>	<b>38.7</b>							± 1.0
	l/l	<b>0.387</b>	<b>0.387</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>250</b>	<b>250</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>10.4</b>	<b>10.6</b>							± 0.5

Ref: TEMP-0387 Rev.47 BACK / VERSO 1300060881