

ABX Minocal



LOT

CX 449

CAL



(Exp.) 2020-09-05

(YYYY - MM - DD)

Rev 1

| PARAMETRES PARAMETERS | UNITES UNITS | WHITEDIFF | | | | | | | | TOLERANCES TOLERANCE |
|--------------------------|---|--------------|--------------------|--|--|--|--|--|--|-------------------------|
| | | YUMIZEN | YUMIZEN | | | | | | | |
| | | H550 | H500 OT H500 CT | | | | | | | |
| GB WBC | 10 ³ /mm ³ ;10 ⁹ /l | 8.70 | 8.60 | | | | | | | ± 0.20 |
| GR RBC | 10 ⁶ /mm ³ ;10 ¹² /l | 4.52 | 4.52 | | | | | | | ± 0.06 |
| HB HGB | g/dl | 13.3 | 13.2 | | | | | | | ± 0.2 |
| | g/l | 133 | 132 | | | | | | | ± 2 |
| | mmol/l | 8.26 | 8.20 | | | | | | | ± 0.12 |
| HT HCT | % | 38.6 | 38.4 | | | | | | | ± 1.0 |
| | l/l | 0.386 | 0.384 | | | | | | | ± 0.010 |
| PLAQ. PLTS | 10 ³ /mm ³ ;10 ⁹ /l | 245 | 240 | | | | | | | ± 10 |
| VPM MPV | µm ³ ; fl | 10.2 | 10.7 | | | | | | | ± 0.5 |

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060881