

ABX Minocal

LOT		CX 448		CAL		(Exp.) 2020-08-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H550	H500 OT H500 CT								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.90	8.90								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40	4.40								± 0.06
HB HGB	g/dl	13.2	13.2								± 0.2
	g/l	132	132								± 2
	mmol/l	8.20	8.20								± 0.12
HT HCT	%	38.1	38.1								± 1.0
	l/l	0.381	0.381								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	258	253								± 10
VPM MPV	µm ³ ; fl	10.3	11.0								± 0.5

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060881