

LOT		CX 426		CAL		(Exp.) 2018-10-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS		UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
			YUMIZEN	YUMIZEN							
			H550	H500 OT H500 CT							
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.1							± 0.2
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.52							± 0.06
HB	HGB	g/dl	13.0	13.0							± 0.2
		g/l	130	130							± 2
		mmol/l	8.07	8.07							± 0.12
HT	HCT	%	37.6	37.7							± 1.0
		l/l	0.376	0.377							± 0.010
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	251	249							± 10
VPM	MPV	µm ³ ; fl	11.0	10.1							± 0.5

Ref: TEMP-0387 Rev.45 BACK / VERSO 1300032436-A