

# ABX Minocal



**LOT**

**CX 426**

**CAL**



**(Exp.) 2018-10-05**

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.9</b>	<b>9.1</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.50</b>	<b>4.52</b>							± 0.06
HB HGB	g/dl	<b>13.0</b>	<b>13.0</b>							± 0.2
	g/l	<b>130</b>	<b>130</b>							± 2
	mmol/l	<b>8.07</b>	<b>8.07</b>							± 0.12
HT HCT	%	<b>37.6</b>	<b>37.7</b>							± 1.0
	l/l	<b>0.376</b>	<b>0.377</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>251</b>	<b>249</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>11.0</b>	<b>10.1</b>							± 0.5

Ref: TEMP-0387 Rev.45 BACK / VERSO 1300032436-A