

ABX Minocal



LOT

CX 457
Rev 1

CAL

(Exp.) **2021-05-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.75	8.75							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.60	4.60							± 0.06
HB HGB	g/dl	13.4	13.4							± 0.2
	g/l	134	134							± 2
	mmol/l	8.32	8.32							± 0.12
HT HCT	%	38.9	38.6							± 1.0
	l/l	0.389	0.386							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	260	253							± 10
VPM MPV	µm ³ ; fl	11.3	11.8							± 0.5

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060881