

ABX Minocal



LOT

CX 437
Rev 1

CAL

(Exp.) **2019-09-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49	4.49							± 0.06
HB HGB	g/dl	13.2	13.2							± 0.2
	g/l	132	132							± 2
	mmol/l	8.20	8.20							± 0.12
HT HCT	%	38.2	37.9							± 1.0
	l/l	0.382	0.379							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	255	258							± 10
VPM MPV	µm ³ ; fl	10.2	10.1							± 0.5

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060881