

# ABX Minocal



**LOT**

**CX 075**

**CAL**



**(Exp.) 2015-08-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN								
		H500								
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.0</b>								± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.51</b>								± 0.06
HB HGB	g/dl	<b>13.7</b>								± 0.2
	g/l	<b>137</b>								± 2
	mmol/l	<b>8.51</b>								± 0.12
HT HCT	%	<b>38.6</b>								± 1.0
	l/l	<b>0.386</b>								± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>233</b>								± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>9.7</b>								± 0.5

Ref: TEMP-0387 Rev.40 BACK / VERSO 1219930095-A