

ABX Minocal



LOT

CX 075

CAL



(Exp.) 2015-08-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE	
		YUMIZEN									
		H500									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0									± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51									± 0.06
HB HGB	g/dl	13.7									± 0.2
	g/l	137									± 2
	mmol/l	8.51									± 0.12
HT HCT	%	38.6									± 1.0
	l/l	0.386									± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	233									± 10
VPM MPV	µm ³ ; fl	9.7									± 0.5

Ref: TEMP-0387 Rev.40 BACK / VERSO 1219930095-A