

ABX Minocal



LOT

CX 418

CAL



(Exp.) 2018-02-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.1							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.54							± 0.06
HB HGB	g/dl	13.4	13.4							± 0.2
	g/l	134	134							± 2
	mmol/l	8.32	8.32							± 0.12
HT HCT	%	39.0	38.8							± 1.0
	l/l	0.390	0.388							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	272	266							± 10
VPM MPV	µm ³ ; fl	10.9	9.9							± 0.5

Ref: TEMP-0387 Rev.44 BACK / VERSO 1300032436-A