

ABX Minocal



LOT

CX 408

CAL



(Exp.) 2017-04-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.1							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.58							± 0.06
HB HGB	g/dl		13.1							± 0.2
	g/l		131							± 2
	mmol/l		8.14							± 0.12
HT HCT	%		37.9							± 1.0
	l/l		0.379							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		267							± 10
VPM MPV	µm ³ ; fl		10.6							± 0.5

Ref: TEMP-0387 Rev.44 BACK / VERSO 1300032436-A