

ABX Minocal



LOT

CX 417
Rev 1

CAL

(Exp.) **2018-01-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.2							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.72	4.74							± 0.06
HB HGB	g/dl	13.2	13.3							± 0.2
	g/l	132	133							± 2
	mmol/l	8.20	8.26							± 0.12
HT HCT	%	39.9	38.5							± 1.0
	l/l	0.399	0.385							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	268	260							± 10
VPM MPV	μm ³ ; fl	10.7	10.3							± 0.5

Ref: TEMP-0387 Rev.44 BACK / VERSO 1300032436-A