

ABX Minocal



LOT

CX 407

CAL



(Exp.) 2017-03-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H500 OT v1 only	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3	9.2							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.62	4.63							± 0.06
HB HGB	g/dl	13.3	13.3							± 0.2
	g/l	133	133							± 2
	mmol/l	8.26	8.26							± 0.12
HT HCT	%	38.1	38.4							± 1.0
	l/l	0.381	0.384							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	251	255							± 10
VPM MPV	µm ³ ; fl	10.4	10.3							± 0.5

Ref: TEMP-0387 Rev.43 BACK / VERSO 1300027265-A