

# ABX Minocal



**LOT**

**CX 406**

**CAL**



**(Exp.) 2017-02-05**

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H500 OT v1 only	H500 OT H500 CT							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.0</b>	<b>9.0</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.65</b>	<b>4.68</b>							± 0.06
HB HGB	g/dl	<b>13.3</b>	<b>13.3</b>							± 0.2
	g/l	<b>133</b>	<b>133</b>							± 2
	mmol/l	<b>8.26</b>	<b>8.26</b>							± 0.12
HT HCT	%	<b>37.9</b>	<b>38.4</b>							± 1.0
	l/l	<b>0.379</b>	<b>0.384</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>242</b>	<b>250</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>10.7</b>	<b>10.7</b>							± 0.5

Ref: TEMP-0387 Rev.43 BACK / VERSO 1300027265-A