

# ABX Minocal

**LOT**

**CX 036**

**CAL**



**(Exp.) 2016-04-05**  
(YYYY - MM - DD)

Rev **2**

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN		YUMIZEN		YUMIZEN		YUMIZEN		
		H500 OT v1 only	H500 OT H500 CT							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8,8</b>	<b>9,0</b>							± 0,2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4,47</b>	<b>4,46</b>							± 0,06
HB HGB	g/dl	<b>13,2</b>	<b>13,2</b>							± 0,2
	g/l	<b>132</b>	<b>132</b>							± 2
	mmol/l	<b>8,20</b>	<b>8,20</b>							± 0,12
HT HCT	%	<b>37,5</b>	<b>38,2</b>							± 1,0
	l/l	<b>0,375</b>	<b>0,382</b>							± 0,010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>271</b>	<b>251</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>10,0</b>	<b>10,0</b>							± 0,5

Ref: TEMP-0387 Rev.42 BACK / VERSO 1300019741-A