

# ABX Minocal



**LOT**

**CX 016**

**CAL**



**(Exp.) 2016-02-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN								
		H500								
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.0</b>								± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.60</b>								± 0.06
HB HGB	g/dl	<b>13.5</b>								± 0.2
	g/l	<b>135</b>								± 2
	mmol/l	<b>8.38</b>								± 0.12
HT HCT	%	<b>39.1</b>								± 1.0
	l/l	<b>0.391</b>								± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>260</b>								± 10
VPM MPV	μm <sup>3</sup> ; fl	<b>9.6</b>								± 0.5

Ref: TEMP-0387 Rev.40 BACK / VERSO 1219930095-A