

LOT		CX 125		CAL		(Exp.) 2016-01-05 (YYYY - MM - DD)					
Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN									
		H500									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1									± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.41									± 0.06
HB HGB	g/dl	13.2									± 0.2
	g/l	132									± 2
	mmol/l	8.20									± 0.12
HT HCT	%	38.2									± 1.0
	l/l	0.382									± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	229									± 10
VPM MPV	µm ³ ; fl	9.9									± 0.5

Ref: TEMP-0387 Rev.40 BACK / VERSO 1219930095-A