

# ABX Minocal



**LOT**

**CX 125**

**CAL**



**(Exp.) 2016-01-05**

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN								
		H500								
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.1</b>								± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.41</b>								± 0.06
HB HGB	g/dl	<b>13.2</b>								± 0.2
	g/l	<b>132</b>								± 2
	mmol/l	<b>8.20</b>								± 0.12
HT HCT	%	<b>38.2</b>								± 1.0
	l/l	<b>0.382</b>								± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>229</b>								± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>9.9</b>								± 0.5

Ref: TEMP-0387 Rev.40 BACK / VERSO 1219930095-A