

ABX Minocal



LOT

CX 115

CAL



(Exp.) 2015-12-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN								
		H500								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56								± 0.06
HB HGB	g/dl	13.3								± 0.2
	g/l	133								± 2
	mmol/l	8.26								± 0.12
HT HCT	%	38.6								± 1.0
	l/l	0.386								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	251								± 10
VPM MPV	µm ³ ; fl	10.1								± 0.5

Ref: TEMP-0387 Rev.40 BACK / VERSO 1219930095-A