

ABX Minocal

LOT		CX 429		CAL		(Exp.) 2019-01-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS		UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
			YUMIZEN	YUMIZEN							
			H550	H500 OT H500 CT							
GB	WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.0							± 0.2
GR	RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.63	4.62							± 0.06
HB	HGB	g/dl	13.6	13.6							± 0.2
		g/l	136	136							± 2
		mmol/l	8.45	8.45							± 0.12
HT	HCT	%	39.1	39.0							± 1.0
		l/l	0.391	0.390							± 0.010
PLAQ.	PLTS	10 ³ /mm ³ ;10 ⁹ /l	252	250							± 10
VPM	MPV	µm ³ ; fl	10.2	9.2							± 0.5

 1300032436-A
 BACK / VERSO
 Ref. TEMP-0387 Rev.45