

# ABX Minocal



**LOT**

**CX 429**

**CAL**



**(Exp.) 2019-01-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.9</b>	<b>9.0</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.63</b>	<b>4.62</b>							± 0.06
HB HGB	g/dl	<b>13.6</b>	<b>13.6</b>							± 0.2
	g/l	<b>136</b>	<b>136</b>							± 2
	mmol/l	<b>8.45</b>	<b>8.45</b>							± 0.12
HT HCT	%	<b>39.1</b>	<b>39.0</b>							± 1.0
	l/l	<b>0.391</b>	<b>0.390</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>252</b>	<b>250</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>10.2</b>	<b>9.2</b>							± 0.5

Ref: TEMP-0387 Rev.45 BACK / VERSO 1300032436-A