

ABX Minocal



LOT

CX 431
Rev 1

CAL

(Exp.) **2019-03-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.60	8.60							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.51							± 0.06
HB HGB	g/dl	12.9	12.9							± 0.2
	g/l	129	129							± 2
	mmol/l	8.01	8.01							± 0.12
HT HCT	%	37.3	37.0							± 1.0
	l/l	0.373	0.370							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	262	257							± 10
VPM MPV	µm ³ ; fl	11.5	10.6							± 0.5

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060881