

ABX Minocal



LOT

CX 442

CAL



(Exp.) 2020-02-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.00	9.00							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.62	4.62							± 0.06
HB HGB	g/dl	13.7	13.7							± 0.2
	g/l	137	137							± 2
	mmol/l	8.51	8.51							± 0.12
HT HCT	%	39.5	39.3							± 1.0
	l/l	0.395	0.393							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	260	255							± 10
VPM MPV	µm ³ ; fl	9.8	9.8							± 0.5

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060881