

ABX Minocal



LOT

CX 422
Rev 1

CAL

(Exp.) **2018-06-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.0							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.46							± 0.06
HB HGB	g/dl	12.5	12.5							± 0.2
	g/l	125	125							± 2
	mmol/l	7.76	7.76							± 0.12
HT HCT	%	36.4	36.0							± 1.0
	l/l	0.364	0.360							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	268	263							± 10
VPM MPV	µm ³ ; fl	10.9	9.8							± 0.5

Ref: TEMP-0387 Rev.45 BACK / VERSO 1300032436-A