

ABX Minocal



LOT

CX 412
Rev 1

CAL

(Exp.) **2017-08-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.1							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.57							± 0.06
HB HGB	g/dl	13.5	13.4							± 0.2
	g/l	135	134							± 2
	mmol/l	8.38	8.32							± 0.12
HT HCT	%	39.8	38.6							± 1.0
	l/l	0.398	0.386							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	247	245							± 10
VPM MPV	µm ³ ; fl	10.2	9.8							± 0.5

Ref: TEMP-0387 Rev.44 BACK / VERSO 1300032436-A