

ABX Minocal



LOT

CX 412
Rev 1

CAL

(Exp.) **2017-08-05**
(YYYY - MM - DD)

| PARAMETRES PARAMETERS | UNITES UNITS | WHITEDIFF | | | | | | | | TOLERANCES TOLERANCE |
|--------------------------|---|--------------|--------------------|--|--|--|--|--|--|-------------------------|
| | | YUMIZEN | YUMIZEN | | | | | | | |
| | | H550 | H500 OT H500 CT | | | | | | | |
| GB WBC | 10 ³ /mm ³ ;10 ⁹ /l | 9.0 | 9.1 | | | | | | | ± 0.2 |
| GR RBC | 10 ⁶ /mm ³ ;10 ¹² /l | 4.58 | 4.57 | | | | | | | ± 0.06 |
| HB HGB | g/dl | 13.5 | 13.4 | | | | | | | ± 0.2 |
| | g/l | 135 | 134 | | | | | | | ± 2 |
| | mmol/l | 8.38 | 8.32 | | | | | | | ± 0.12 |
| HT HCT | % | 39.8 | 38.6 | | | | | | | ± 1.0 |
| | l/l | 0.398 | 0.386 | | | | | | | ± 0.010 |
| PLAQ. PLTS | 10 ³ /mm ³ ;10 ⁹ /l | 247 | 245 | | | | | | | ± 10 |
| VPM MPV | µm ³ ; fl | 10.2 | 9.8 | | | | | | | ± 0.5 |

Ref: TEMP-0387 Rev.44 BACK / VERSO 1300032436-A