

ABX Minocal


LOT
CX 402
CAL

(Exp.) 2016-10-05

(YYYY - MM - DD)

 Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H500 OT v1 only	H500 OT H500 CT								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.2								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.47								± 0.06
HB HGB	g/dl	13.2	13.4								± 0.2
	g/l	132	134								± 2
	mmol/l	8.20	8.32								± 0.12
HT HCT	%	37.2	38.4								± 1.0
	l/l	0.372	0.384								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	252	243								± 10
VPM MPV	µm ³ ; fl	10.5	11.0								± 0.5

Ref: TEMP-0387 Rev.42 BACK / VERSO 1300019741-A