

ABX Minocal



LOT

CX 451

CAL



(Exp.) 2020-11-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.00	9.00							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.55							± 0.06
HB HGB	g/dl	13.6	13.5							± 0.2
	g/l	136	135							± 2
	mmol/l	8.45	8.38							± 0.12
HT HCT	%	39.2	39.2							± 1.0
	l/l	0.392	0.392							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	258	252							± 10
VPM MPV	µm ³ ; fl	10.9	11.5							± 0.5

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060881