

ABX Minocal

LOT		CX 441		CAL		(Exp.) 2020-01-05 (YYYY - MM - DD)					
Rev 1											
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H550	H500 OT H500 CT								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.85	8.90								± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.58								± 0.06
HB HGB	g/dl	13.5	13.5								± 0.2
	g/l	135	135								± 2
	mmol/l	8.38	8.38								± 0.12
HT HCT	%	38.8	38.8								± 1.0
	l/l	0.388	0.388								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	260	260								± 10
VPM MPV	µm ³ ; fl	11.0	11.0								± 0.5

Ref: TEMP-0387 Rev.46

1300060881

BACK / VERSO

 1300060881
 BACK / VERSO
 Ref: TEMP-0387 Rev.46