

ABX Minocal



LOT

CX 421
Rev 1

CAL

(Exp.) **2018-05-05**
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.7	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.88	4.86							± 0.06
HB HGB	g/dl	13.5	13.5							± 0.2
	g/l	135	135							± 2
	mmol/l	8.38	8.38							± 0.12
HT HCT	%	39.3	39.1							± 1.0
	l/l	0.393	0.391							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	265	250							± 10
VPM MPV	µm ³ ; fl	11.3	10.2							± 0.5

Ref: TEMP-0387 Rev.45 BACK / VERSO 1300032436-A