

# ABX Minocal



**LOT**

**CX 411**  
Rev 1

**CAL**

(Exp.) **2017-07-05**  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.9</b>	<b>9.0</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.65</b>	<b>4.67</b>							± 0.06
HB HGB	g/dl	<b>13.6</b>	<b>13.4</b>							± 0.2
	g/l	<b>136</b>	<b>134</b>							± 2
	mmol/l	<b>8.45</b>	<b>8.32</b>							± 0.12
HT HCT	%	<b>39.4</b>	<b>38.8</b>							± 1.0
	l/l	<b>0.394</b>	<b>0.388</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>265</b>	<b>270</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>10.5</b>	<b>10.0</b>							± 0.5

Ref: TEMP-0387 Rev.44 BACK / VERSO 1300032436-A