

ABX Minocal



LOT

CX 401
Rev 1

CAL

(Exp.) 2016-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H500 OT v1 only	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.57	4.59							± 0.06
HB HGB	g/dl	13.3	13.4							± 0.2
	g/l	133	134							± 2
	mmol/l	8.26	8.32							± 0.12
HT HCT	%	37.3	38.4							± 1.0
	l/l	0.373	0.384							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	265	259							± 10
VPM MPV	μm ³ ; fl	9.0	9.6							± 0.5

Ref: TEMP-0387 Rev.42 BACK / VERSO 1300019741-A