

# ABX Minocal



**LOT**

**CX 450**

**CAL**



**(Exp.) 2020-10-05**

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.70</b>	<b>8.70</b>							± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.54</b>	<b>4.54</b>							± 0.06
HB HGB	g/dl	<b>13.4</b>	<b>13.4</b>							± 0.2
	g/l	<b>134</b>	<b>134</b>							± 2
	mmol/l	<b>8.32</b>	<b>8.32</b>							± 0.12
HT HCT	%	<b>39.0</b>	<b>39.0</b>							± 1.0
	l/l	<b>0.390</b>	<b>0.390</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>260</b>	<b>255</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>11.2</b>	<b>11.8</b>							± 0.5

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060881