

ABX Minocal



LOT

CX 430

CAL



(Exp.) 2019-02-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.63	4.61							± 0.06
HB HGB	g/dl	13.2	13.2							± 0.2
	g/l	132	132							± 2
	mmol/l	8.20	8.20							± 0.12
HT HCT	%	38.4	37.9							± 1.0
	l/l	0.384	0.379							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	268	264							± 10
VPM MPV	µm ³ ; fl	11.4	10.5							± 0.5

Ref: TEMP-0387 Rev.45 BACK / VERSO 1300032436-A