

ABX Minocal



LOT

CX 420

CAL



(Exp.) 2018-04-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.7	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.62	4.62							± 0.06
HB HGB	g/dl	13.6	13.5							± 0.2
	g/l	136	135							± 2
	mmol/l	8.45	8.38							± 0.12
HT HCT	%	38.9	38.8							± 1.0
	l/l	0.389	0.388							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	263	254							± 10
VPM MPV	μm ³ ; fl	10.8	9.8							± 0.5

Ref: TEMP-0387 Rev.45 BACK / VERSO 1300032436-A