

ABX Minocal



LOT		CX 400		CAL		(Exp.) 2016-08-05 (YYYY - MM - DD)					
PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF									TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN								
		H500 OT v1 only	H500 OT H500 CT								
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	8.9								± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.72	4.68								± 0.06
HB HGB	g/dl	13.4	13.4								± 0.2
	g/l	134	134								± 2
	mmol/l	8.32	8.32								± 0.12
HT HCT	%	38.5	38.8								± 1.0
	l/l	0.385	0.388								± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	270	253								± 10
VPM MPV	µm ³ ; fl	9.8	10.1								± 0.5

Ref: TEMP-0387 Rev.42 BACK / VERSO 1300019741-A