

ABX Minocal



LOT

CX 400
Rev 1

CAL

(Exp.) 2016-08-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H500 OT v1 only	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.72	4.68							± 0.06
HB HGB	g/dl	13.4	13.4							± 0.2
	g/l	134	134							± 2
	mmol/l	8.32	8.32							± 0.12
HT HCT	%	38.5	38.8							± 1.0
	l/l	0.385	0.388							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	270	253							± 10
VPM MPV	µm ³ ; fl	9.8	10.1							± 0.5

Ref: TEMP-0387 Rev.42 BACK / VERSO 1300019741-A