



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT MX 054 L
 Rev 1

CONTROL L



2014-08-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA		PENTRA
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		1.9	1.9		1.9					± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		2.30	2.34		2.30					± 0.15
HB HGB	g/dl		5.8	6.1		5.6					± 0.4
	g/l		58	61		56					± 4
	mmol/l		3.59	3.79		3.47					± 0.25
HT HCT	%		15.9	16.4		16.1					± 2.0
	l/l		0.159	0.164		0.161					± 0.020
VGM MCV	µm ³ ; fl		69	70		70					± 5
TGMH MCH	pg		25.1	26.1		24.3					± 2.0
	fmol		1.56	1.62		1.51					± 0.12
CCMH MCHC	g/dl		36.4	37.2		34.7					± 3.0
	g/l		364	372		347					± 30
	mmol/l		22.63	23.13		21.52					± 1.86
IDR RDW	%		13.8	13.2		13.8					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		77	78		75					± 20
VPM MPV	µm ³ ; fl		8.7	8.4		8.1					± 2.0
LYMPHO	#		1.24	1.2		1.19					± 0.40
	%		65.1	61.0		62.5					± 7.0
MONO	#		0.19	0.2		0.18					± 0.17
	%		10.0	9.5		9.5					± 6.0
GRANULO	#		0.47	0.6		0.53					± 0.40
	%		25.0	29.5		28.0					± 7.0

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 054 N
 Rev 1

CONTROL N

(Exp.)

2014-08-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA		PENTRA
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		7.5	7.8		7.4					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		4.50	4.45		4.50					± 0.18
HB HGB	g/dl		13.3	13.1		13.0					± 0.5
	g/l		133	131		130					± 5
	mmol/l		8.26	8.14		8.09					± 0.31
HT HCT	%		36.9	36.0		36.9					± 2.5
	l/l		0.369	0.360		0.369					± 0.025
VGM MCV	µm ³ ; fl		82	81		82					± 5
TGMH MCH	pg		29.6	29.4		29.0					± 2.0
	fmol		1.84	1.83		1.80					± 0.12
CCMH MCHC	g/dl		36.0	36.3		35.3					± 3.0
	g/l		360	363		353					± 30
	mmol/l		22.38	22.57		21.94					± 1.86
IDR RDW	%		13.5	12.8		13.5					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		256	250		255					± 40
VPM MPV	µm ³ ; fl		7.6	7.5		7.2					± 2.0
LYMPHO	#		2.58	2.6		2.44					± 0.40
	%		34.5	33.5		33.2					± 5.0
MONO	#		0.53	0.5		0.51					± 0.30
	%		7.1	7.0		7.0					± 4.0
GRANULO	#		4.37	4.6		4.42					± 0.60
	%		58.3	59.5		60.0					± 6.0

Ref: TEMP-0630 Rev.35 9930086-B



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT

MX 054 H
 Rev 1

CONTROL H

(Exp.)

2014-08-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA		PENTRA
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		20.2	20.6		20.1					± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l		5.48	5.43		5.50					± 0.20
HB HGB	g/dl		17.5	17.4		17.4					± 0.6
	g/l		175	174		174					± 6
	mmol/l		10.87	10.81		10.82					± 0.37
HT HCT	%		48.8	48.3		49.5					± 3.0
	l/l		0.488	0.483		0.495					± 0.030
VGM MCV	µm ³ ; fl		89	89		90					± 5
TGMH MCH	pg		31.9	32.0		31.7					± 2.0
	fmol		1.98	1.99		1.97					± 0.12
CCMH MCHC	g/dl		35.9	36.0		35.2					± 3.0
	g/l		359	360		352					± 30
	mmol/l		22.28	22.36		21.86					± 1.86
IDR RDW	%		13.0	12.1		13.0					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		492	485		495					± 45
VPM MPV	µm ³ ; fl		7.7	7.8		7.3					± 2.0
LYMPHO	#		3.32	3.4		3.10					± 1.00
	%		16.4	16.5		15.5					± 5.0
MONO	#		1.20	1.0		1.00					± 0.80
	%		5.9	5.0		5.0					± 4.0
GRANULO	#		15.64	16.2		15.90					± 1.40
	%		77.4	78.5		79.2					± 7.0