

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 114 L
 Rev 1

CONTROL L

(Exp.)

2015-02-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA		PENTRA
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		CRP 200			60 8P - 60 18P					± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l										± 0.15
HB HGB	g/dl		5.9			5.6					± 0.4
	g/l		59			56					± 4
	mmol/l		3.65			3.47					± 0.25
HT HCT	%		15.7			16.1					± 2.0
	l/l		0.157			0.161					± 0.020
VGM MCV	µm ³ ; fl		69			70					± 4
TGMH MCH	pg		25.9			24.3					± 2.0
	fmol		1.61			1.51					± 0.12
CCMH MCHC	g/dl		37.5			34.7					± 3.0
	g/l		375			347					± 30
	mmol/l		23.31			21.52					± 1.86
IDR RDW	%		13.5			14.0					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		80			85					± 20
VPM MPV	µm ³ ; fl		8.4			7.6					± 2.0
LYMPHO	#		1.32			1.28					± 0.40
	%		66.2			64.0					± 8.0
MONO	#		0.19			0.18					± 0.20
	%		9.5			9.0					± 6.0
GRANULO	#		0.49			0.54					± 0.40
	%		24.5			27.0					± 7.0

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 114 N
 Rev 1

CONTROL N

(Exp.)

2015-02-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA		PENTRA
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		CRP 200			60 8P - 60 18P					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l										± 0.18
HB HGB	g/dl		13.4			13.2					± 0.5
	g/l		134			132					± 5
	mmol/l		8.32			8.22					± 0.31
HT HCT	%		37.2			37.4					± 2.5
	l/l		0.372			0.374					± 0.025
VGM MCV	µm ³ ; fl		81			81					± 4
TGMH MCH	pg		29.2			28.6					± 2.0
	fmol		1.81			1.78					± 0.12
CCMH MCHC	g/dl		36.0			35.4					± 3.0
	g/l		360			354					± 30
	mmol/l		22.38			21.95					± 1.86
IDR RDW	%		13.0			13.5					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		255			255					± 40
VPM MPV	µm ³ ; fl		7.9			7.4					± 2.0
LYMPHO	#		2.69			2.51					± 0.40
	%		35.0			33.7					± 5.0
MONO	#		0.55			0.52					± 0.30
	%		7.1			7.0					± 4.0
GRANULO	#		4.45			4.44					± 0.60
	%		57.8			59.5					± 6.0

9930092-A
 Ref: TEMP-0630 Rev.36

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 114 H
 Rev 1

CONTROL H

(Exp.)

2015-02-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA		PENTRA
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		CRP 200			60 8P - 60 18P					± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l										± 0.20
HB HGB	g/dl		17.7			17.6					± 0.6
	g/l		177			176					± 6
	mmol/l		10.99			10.94					± 0.37
HT HCT	%		49.4			49.8					± 3.0
	l/l		0.494			0.498					± 0.030
VGM MCV	µm ³ ; fl		89			89					± 4
TGMH MCH	pg		31.9			31.5					± 2.0
	fmol		1.98			1.95					± 0.12
CCMH MCHC	g/dl		35.8			35.4					± 3.0
	g/l		358			354					± 30
	mmol/l		22.25			21.96					± 1.86
IDR RDW	%		12.8			13.0					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		525			530					± 45
VPM MPV	µm ³ ; fl		7.7			7.2					± 2.0
LYMPHO	#		3.49			3.31					± 1.00
	%		17.0			16.5					± 5.0
MONO	#		0.89			0.80					± 0.80
	%		4.3			4.0					± 4.0
GRANULO	#		16.15			15.90					± 1.40
	%		78.4			79.2					± 7.0