

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 015 L
 Rev 1

CONTROL L

(Exp.)

2015-04-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA		PENTRA
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		CRP 200			60 8P - 60 18P					± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l										± 0.15
HB HGB	g/dl		6.0			5.7					± 0.4
	g/l		60			57					± 4
	mmol/l		3.71			3.52					± 0.25
HT HCT	%		16.2			16.4					± 2.0
	l/l		0.162			0.164					± 0.020
VGM MCV	µm ³ ; fl		69			70					± 4
TGMH MCH	pg		25.4			24.2					± 2.0
	fmol		1.58			1.51					± 0.12
CCMH MCHC	g/dl		36.9			34.6					± 3.0
	g/l		369			346					± 30
	mmol/l		22.89			21.51					± 1.86
IDR RDW	%		14.0			13.5					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		76			77					± 20
VPM MPV	µm ³ ; fl		8.7			8.0					± 2.0
LYMPHO	#		1.38			1.35					± 0.40
	%		65.7			64.5					± 8.0
MONO	#		0.19			0.20					± 0.20
	%		9.0			9.5					± 6.0
GRANULO	#		0.54			0.55					± 0.40
	%		25.5			26.0					± 7.0

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT MX 015 N
 Rev 1

CONTROL N

(Exp.)

2015-04-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA		PENTRA
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		CRP 200			60 8P - 60 18P					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l										± 0.18
HB HGB	g/dl		7.6			7.5					± 0.5
	g/l		135			133					± 5
	mmol/l		8.38			8.28					± 0.31
HT HCT	%		37.1			37.3					± 2.5
	l/l		0.371			0.373					± 0.025
VGM MCV	µm ³ ; fl		81			81					± 4
TGMH MCH	pg		29.5			29.0					± 2.0
	fmol		1.83			1.80					± 0.12
CCMH MCHC	g/dl		36.4			35.8					± 3.0
	g/l		364			358					± 30
	mmol/l		22.60			22.21					± 1.86
IDR RDW	%		13.0			13.3					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		255			260					± 40
VPM MPV	µm ³ ; fl		8.3			7.6					± 2.0
LYMPHO	#		2.62			2.51					± 0.40
	%		34.5			33.7					± 5.0
MONO	#		0.51			0.49					± 0.30
	%		6.7			6.5					± 4.0
GRANULO	#		4.47			4.48					± 0.60
	%		58.8			60.0					± 6.0

9930092-A
 Ref: TEMP-0630 Rev.36



HORIBA ABX SAS

HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17



LOT MX 015 H
 Rev 1

CONTROL H

(Exp.)

2015-04-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	ARGOS - HELIOS	PENTRA	PENTRA		PENTRA
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l		CRP 200			60 8P - 60 18P					± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l										± 0.20
HB HGB	g/dl		17.7			17.7					± 0.6
	g/l		177			177					± 6
	mmol/l		10.99			11.00					± 0.37
HT HCT	%		49.5			50.7					± 3.0
	l/l		0.495			0.507					± 0.030
VGM MCV	µm ³ ; fl		88			89					± 4
TGMH MCH	pg		31.4			31.1					± 2.0
	fmol		1.95			1.93					± 0.12
CCMH MCHC	g/dl		35.7			34.9					± 3.0
	g/l		357			349					± 30
	mmol/l		22.19			21.69					± 1.86
IDR RDW	%		12.9			13.0					± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l		505			520					± 45
VPM MPV	µm ³ ; fl		7.8			7.3					± 2.0
LYMPHO	#		3.04			2.88					± 1.00
	%		14.8			14.4					± 5.0
MONO	#		1.33			1.10					± 0.80
	%		6.5			5.5					± 4.0
GRANULO	#		16.07			15.92					± 1.40
	%		78.4			79.7					± 7.0