

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT		MX 435 L		CONTROL		L		(Exp.)		2022-08-05			
		Rev 1		LyseVet									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS								TOLERANCES TOLERANCE	
		Sol Vet ABC	Sol Vet ABC Plus ESV 60										
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.7	1.7									± 0.4	
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	2.26	2.24									± 0.15	
HB HGB	g/dl	5.8	5.7									± 0.4	
	g/l	58	57									± 4	
	mmol/l	3.60	3.54									± 0.25	
HT HCT	%	16.0	15.9									± 2.0	
	l/l	0.160	0.159									± 0.020	
VGM MCV	µm ³ ; fl	71	71									± 4	
TGMH MCH	pg	25.7	25.4									± 2.0	
	fmol	1.59	1.58									± 0.12	
	g/dl	36.1	35.8									± 3.0	
CCMH MCHC	g/l	361	358									± 30	
	mmol/l	22.45	22.26									± 1.86	
	%	N/A	15.0									± 3.0	
IDR RDW	%	N/A	15.0									± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	71	75									± 20	
VPM MPV	µm ³ ; fl	N/A	9.1									± 2.0	
LYMPHO	#	N/A	0.93									± 0.40	
	%	N/A	54.5									± 8.0	
MONO	#	N/A	0.21									± 0.20	
	%	N/A	12.5									± 6.0	
GRANULO	#	N/A	0.56									± 0.40	
	%	N/A	33.0									± 7.0	

LOT		MX 435 N		CONTROL		N		(Exp.)		2022-08-05			
		Rev 1		LyseVet									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS								TOLERANCES TOLERANCE	
		Sol Vet ABC	Sol Vet ABC Plus ESV 60										
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.0	7.3									± 0.8	
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	4.63	4.63									± 0.18	
HB HGB	g/dl	13.4	13.4									± 0.5	
	g/l	134	134									± 5	
	mmol/l	8.32	8.32									± 0.31	
HT HCT	%	37.0	37.5									± 2.5	
	l/l	0.370	0.375									± 0.025	
VGM MCV	µm ³ ; fl	80	81									± 4	
TGMH MCH	pg	28.9	28.9									± 2.0	
	fmol	1.80	1.80									± 0.12	
	g/dl	36.2	35.7									± 3.0	
CCMH MCHC	g/l	362	357									± 30	
	mmol/l	22.47	22.19									± 1.86	
	%	N/A	14.5									± 3.0	
IDR RDW	%	N/A	14.5									± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	245	255									± 40	
VPM MPV	µm ³ ; fl	N/A	8.8									± 2.0	
LYMPHO	#	N/A	2.04									± 0.40	
	%	N/A	28.0									± 5.0	
MONO	#	N/A	0.62									± 0.30	
	%	N/A	8.5									± 4.0	
GRANULO	#	N/A	4.64									± 0.60	
	%	N/A	63.5									± 6.0	

Ref. TEMP-050 Rev.41 BACK / VERSO 1300060626

LOT		MX 435 H		CONTROL		H		(Exp.)		2022-08-05			
		Rev 1		LyseVet									
PARAMETRES PARAMETERS	UNITES UNITS	MICROS		MICROS								TOLERANCES TOLERANCE	
		Sol Vet ABC	Sol Vet ABC Plus ESV 60										
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	18.8	19.0									± 1.6	
GR RBC	10 ¹² /mm ³ ; 10 ¹² /l	5.65	5.65									± 0.20	
HB HGB	g/dl	18.1	18.0									± 0.6	
	g/l	181	180									± 6	
	mmol/l	11.24	11.18									± 0.37	
HT HCT	%	50.3	50.3									± 3.0	
	l/l	0.503	0.503									± 0.030	
VGM MCV	µm ³ ; fl	89	89									± 4	
TGMH MCH	pg	32.0	31.9									± 2.0	
	fmol	1.99	1.98									± 0.12	
	g/dl	36.0	35.8									± 3.0	
CCMH MCHC	g/l	360	358									± 30	
	mmol/l	22.35	22.23									± 1.86	
	%	N/A	13.5									± 3.0	
IDR RDW	%	N/A	13.5									± 3.0	
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	500	525									± 45	
VPM MPV	µm ³ ; fl	N/A	8.5									± 2.0	
LYMPHO	#	N/A	2.09									± 1.00	
	%	N/A	11.0									± 5.0	
MONO	#	N/A	0.95									± 0.80	
	%	N/A	5.0									± 4.0	
GRANULO	#	N/A	15.96									± 1.40	
	%	N/A	84.0									± 7.0	