

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

MX 409 L
Rev 1

CONTROL L

 (Exp.)

2018-04-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	1.9							± 0.4
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	2.28	2.27							± 0.15
HB HGB	g/dl	5.9	5.9							± 0.4
	g/l	59	59							± 4
	mmol/l	3.66	3.66							± 0.25
HT HCT	%	16.2	15.9							± 2.0
	l/l	0.162	0.159							± 0.020
VGM MCV	µm ³ ; fl	71	70							± 4
TGMH MCH	pg	25.9	26.0							± 2.0
	fmol	1.61	1.61							± 0.12
CCMH MCHC	g/dl	36.4	37.1							± 3.0
	g/l	364	371							± 30
	mmol/l	22.63	23.06							± 1.86
IDR RDW	%	N/A	14.0							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	80	83							± 20
VPM MPV	µm ³ ; fl	N/A	8.1							± 2.0
LYMPHO	#	N/A	1.03							± 0.40
	%	N/A	54.0							± 8.0
MONO	#	N/A	0.20							± 0.20
	%	N/A	10.0							± 6.0
GRANULO	#	N/A	0.68							± 0.40
	%	N/A	36.0							± 7.0

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT

MX 409 N
 Rev 1

CONTROL N

(Exp.)

2018-04-05
 (YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.1	7.3							± 0.8
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	4.72	4.72							± 0.18
HB HGB	g/dl	13.6	13.6							± 0.5
	g/l	136	136							± 5
	mmol/l	8.45	8.45							± 0.31
HT HCT	%	37.8	37.3							± 2.5
	l/l	0.378	0.373							± 0.025
VGM MCV	µm ³ ; fl	80	79							± 4
TGMH MCH	pg	28.8	28.8							± 2.0
	fmol	1.79	1.79							± 0.12
CCMH MCHC	g/dl	36.0	36.5							± 3.0
	g/l	360	365							± 30
	mmol/l	22.37	22.65							± 1.86
IDR RDW	%	N/A	14.0							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	245	249							± 40
VPM MPV	µm ³ ; fl	N/A	8.9							± 2.0
LYMPHO	#	N/A	2.12							± 0.40
	%	N/A	29.0							± 5.0
MONO	#	N/A	0.58							± 0.30
	%	N/A	8.0							± 4.0
GRANULO	#	N/A	4.60							± 0.60
	%	N/A	63.0							± 6.0

1300010912-A
 Ref: TEMP-0830 Rev.38

HORIBA ABX SAS
 HORIBA ABX SAS
 Parc Euromédecine
 Rue du Caducée - BP 7290
 34184 MONTPELLIER Cedex 4
 FRANCE
 Tél.: 33 (0) 4 67 14 15 16
 Fax.: 33 (0) 4 67 14 15 17

LOT

MX 409 H
 Rev 1

CONTROL H

(Exp.)

2018-04-05
 (YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.5	20.0							± 1.6
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	5.70	5.74							± 0.20
HB HGB	g/dl	18.0	17.9							± 0.6
	g/l	180	179							± 6
	mmol/l	11.18	11.12							± 0.37
HT HCT	%	50.2	49.9							± 3.0
	l/l	0.502	0.499							± 0.030
VGM MCV	µm ³ ; fl	88	87							± 4
TGMH MCH	pg	31.6	31.2							± 2.0
	fmol	1.96	1.94							± 0.12
CCMH MCHC	g/dl	35.9	35.8							± 3.0
	g/l	359	358							± 30
	mmol/l	22.28	22.26							± 1.86
IDR RDW	%	N/A	13.5							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	495	488							± 45
VPM MPV	µm ³ ; fl	N/A	8.1							± 2.0
LYMPHO	#	N/A	2.30							± 1.00
	%	N/A	11.5							± 5.0
MONO	#	N/A	1.00							± 0.80
	%	N/A	5.0							± 4.0
GRANULO	#	N/A	16.70							± 1.40
	%	N/A	83.5							± 7.0