

ABX Minotrol 16



LOT MX 112 L

Rev 1

CONTROL L

(Exp.) 2013-02-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK						TOLERANCES TOLERANCE
		MICROS		MICROS				
		Scil Vet ABC	Scil Vet ABC Plus ESV 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.7	1.8					± 0.4
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	2.32	2.31					± 0.15
HB HGB	g/dl	5.7	5.8					± 0.4
	g/l	57	58					± 4
	mmol/l	3.54	3.60					± 0.25
HT HCT	%	16.0	16.2					± 2.0
	l/l	0.160	0.162					± 0.020
VGM MCV	µm ³ ;fl	69	70					± 4
TGMH MCH	pg	24.6	25.1					± 2.0
	f mol	1.53	1.56					± 0.12
CCMH MCHC	g/dl	35.6	35.9					± 3.0
	g/l	356	359					± 30
	mmol/l	22.11	22.27					± 1.86
IDR RDW	%	N/A	14.0					± 4.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	64	72					± 20
VPM MPV	µm ³ ;fl	N/A	8.0					± 2.0
LYMPHO	#	N/A	1.0					± 0.40
	%	N/A	54.0					± 9.0
MONO	#	N/A	0.2					± 0.20
	%	N/A	12.0					± 7.0
GRANULO	#	N/A	0.6					± 0.40
	%	N/A	34.0					± 8.0

9930041-E

Ref: TEMP-0830 Rev.17

ABX Minotrol 16



LOT **MX 112 N**
 Rev 1

CONTROL **N**

(Exp.) **2013-02-05**
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK						TOLERANCES TOLERANCE
		MICROS		MICROS				
		Scil Vet ABC	Scil Vet ABC Plus ESV 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.3	7.7					± 0.8
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4.72	4.70					± 0.18
HB HGB	g/dl	13.4	13.4					± 0.5
	g/l	134	134					± 5
	mmol/l	8.32	8.32					± 0.31
HT HCT	%	37.8	38.1					± 2.5
	l/l	0.378	0.381					± 0.025
VGM MCV	µm ³ ; fl	80	81					± 4
TGMH MCH	pg	28.4	28.5					± 2.0
	f mol	1.76	1.77					± 0.12
CCMH MCHC	g/dl	35.5	35.2					± 3.0
	g/l	355	352					± 30
	mmol/l	22.04	21.86					± 1.86
IDR RDW	%	N/A	13.5					± 4.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	258	260					± 40
VPM MPV	µm ³ ; fl	N/A	7.5					± 2.0
LYMPHO	#	N/A	2.0					± 0.50
	%	N/A	26.5					± 6.0
MONO	#	N/A	0.7					± 0.40
	%	N/A	9.0					± 5.0
GRANULO	#	N/A	5.0					± 0.60
	%	N/A	64.5					± 7.0

9930042-E

Ref: TEMP-0830 Rev.17

ABX Minotrol 16



LOT **MX 112 H**
 Rev 1

CONTROL **H**

(Exp.) **2013-02-05**
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	19.5	20.2							± 1.6
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	5.77	5.74							± 0.20
HB HGB	g/dl	18.0	18.0							± 0.6
	g/l	180	180							± 6
	mmol/l	11.18	11.18							± 0.37
HT HCT	%	50.8	51.7							± 3.0
	l/l	0.508	0.517							± 0.030
VGM MCV	µm ³ ; fl	88	90							± 4
TGMH MCH	pg	31.2	31.4							± 2.0
	fmol	1.94	1.95							± 0.12
CCMH MCHC	g/dl	35.4	34.8							± 3.0
	g/l	354	348							± 30
	mmol/l	22.01	21.64							± 1.86
IDR RDW	%	N/A	13.0							± 4.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	535	500							± 55
VPM MPV	µm ³ ; fl	N/A	7.7							± 2.0
LYMPHO	#	N/A	2.3							± 1.00
	%	N/A	11.5							± 5.0
MONO	#	N/A	1.0							± 0.80
	%	N/A	5.0							± 4.0
GRANULO	#	N/A	16.9							± 1.40
	%	N/A	83.5							± 7.0

9930043-E

Ref: TEMP-0830 Rev.17