

ABX Minotrol 16



LOT **MX 092 L**

CONTROL **L**

(Exp.) **2012-12-05**
 (YYYY - MM - DD)

Rev 3

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	1.8	2.2							± 0.4
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	2.44	2.40							± 0.15
HB HGB	g/dl	6.1	5.9							± 0.4
	g/l	61	59							± 4
	mmol/l	3.79	3.66							± 0.25
HT HCT	%	17.1	16.6							± 2.0
	l/l	0.171	0.166							± 0.020
VGM MCV	µm ³ ; fl	70	69							± 4
TGMH MCH	pg	25.0	24.6							± 2.0
	f mol	1.55	1.53							± 0.12
CCMH MCHC	g/dl	35.7	35.6							± 3.0
	g/l	357	356							± 30
	mmol/l	22.18	22.13							± 1.86
IDR RDW	%	N/A	13.5							± 4.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	70	72							± 20
VPM MPV	µm ³ ; fl	N/A	7.9							± 2.0
LYMPHO	#	N/A	1.2							± 0.40
	%	N/A	53.0							± 9.0
MONO	#	N/A	0.3							± 0.20
	%	N/A	14.0							± 7.0
GRANULO	#	N/A	0.7							± 0.40
	%	N/A	33.0							± 8.0

9930041-E

Ref: TEMP-0830 Rev.17

ABX Minotrol 16



LOT **MX 092 N**
 Rev 3

CONTROL **N**

(Exp.) **2012-12-05**
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK						TOLERANCES TOLERANCE
		MICROS		MICROS				
		Scil Vet ABC	Scil Vet ABC Plus ESV 60					
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.2	8.2					± 0.8
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	4.73	4.66					± 0.18
HB HGB	g/dl	13.5	13.2					± 0.5
	g/l	135	132					± 5
	mmol/l	8.38	8.20					± 0.31
HT HCT	%	37.8	37.3					± 2.5
	l/l	0.378	0.373					± 0.025
VGM MCV	µm ³ ; fl	80	80					± 4
TGMH MCH	pg	28.5	28.3					± 2.0
	f mol	1.77	1.76					± 0.12
CCMH MCHC	g/dl	35.7	35.4					± 3.0
	g/l	357	354					± 30
	mmol/l	22.16	21.99					± 1.86
IDR RDW	%	N/A	13.4					± 4.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	257	228					± 40
VPM MPV	µm ³ ; fl	N/A	7.4					± 2.0
LYMPHO	#	N/A	2.0					± 0.50
	%	N/A	24.5					± 6.0
MONO	#	N/A	0.7					± 0.40
	%	N/A	8.5					± 5.0
GRANULO	#	N/A	5.5					± 0.60
	%	N/A	67.0					± 7.0

9930042-E

Ref: TEMP-0830 Rev.17

ABX Minotrol 16



LOT MX 092 H

CONTROL H

(Exp.) 2012-12-05
 (YYYY - MM - DD)

Rev 3

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	19.7	22.3							± 1.8
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	5.73	5.70							± 0.22
HB HGB	g/dl	17.8	17.6							± 0.7
	g/l	178	176							± 6
	mmol/l	11.05	10.93							± 0.37
HT HCT	%	50.4	50.7							± 3.0
	l/l	0.504	0.507							± 0.030
VGM MCV	µm ³ ; fl	88	89							± 4
TGMH MCH	pg	31.1	30.9							± 2.0
	f mol	1.93	1.92							± 0.12
CCMH MCHC	g/dl	35.3	34.7							± 3.0
	g/l	353	347							± 30
	mmol/l	21.92	21.54							± 1.86
IDR RDW	%	N/A	12.8							± 4.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	510	430							± 45
VPM MPV	µm ³ ; fl	N/A	7.6							± 2.0
LYMPHO	#	N/A	2.6							± 1.00
	%	N/A	11.5							± 5.0
MONO	#	N/A	1.2							± 0.80
	%	N/A	5.5							± 4.0
GRANULO	#	N/A	18.5							± 1.40
	%	N/A	83.0							± 7.0

9930043-E

Ref: TEMP-0830 Rev.17