

ABX Minotrol 16



LOT **MX 032 L**
 Rev 1

CONTROL **L**

(Exp.) **2012-06-05**
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	1.8	2.0							± 0.4
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	2.35	2.40							± 0.15
HB HGB	g/dl	6.0	6.3							± 0.4
	g/l	60	63							± 4
	mmol/l	3.73	3.89							± 0.25
HT HCT	%	17.2	16.9							± 2.0
	l/l	0.172	0.169							± 0.020
VGM MCV	µm ³ ;fl	73	71							± 4
TGMH MCH	pg	25.5	26.1							± 2.0
	f mol	1.59	1.62							± 0.12
CCMH MCHC	g/dl	35.0	37.0							± 3.0
	g/l	350	370							± 30
	mmol/l	21.72	22.96							± 1.86
IDR RDW	%	N/A	14.6							± 3.0
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	74	96							± 20
VPM MPV	µm ³ ;fl	N/A	7.0							± 2.0
LYMPHO	#	N/A	1.01							± 0.40
	%	N/A	51.0							± 9.0
MONO	#	N/A	0.28							± 0.20
	%	N/A	14.3							± 7.0
GRANULO	#	N/A	0.68							± 0.40
	%	N/A	34.3							± 8.0

Ref: FORM-0630 Rev.16 9930041-E



ABX Minotrol 16



LOT **MX 032 N**
 Rev 1

CONTROL **N**

(Exp.) **2012-06-05**
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	6.8	7.3							± 0.8
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	4.60	4.64							± 0.18
HB HGB	g/dl	13.3	13.1							± 0.5
	g/l	133	131							± 5
	mmol/l	8.26	8.15							± 0.31
HT HCT	%	38.2	37.7							± 2.5
	l/l	0.382	0.377							± 0.025
VGM MCV	µm ³ ; fl	83	81							± 4
TGMH MCH	pg	28.9	28.3							± 2.0
	f mol	1.80	1.76							± 0.12
CCMH MCHC	g/dl	34.8	34.8							± 3.0
	g/l	348	348							± 30
	mmol/l	21.63	21.63							± 1.86
IDR RDW	%	N/A	13.6							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	267	269							± 40
VPM MPV	µm ³ ; fl	N/A	6.9							± 2.0
LYMPHO	#	N/A	2.08							± 0.50
	%	N/A	28.4							± 6.0
MONO	#	N/A	0.64							± 0.40
	%	N/A	8.8							± 5.0
GRANULO	#	N/A	4.55							± 0.60
	%	N/A	62.1							± 7.0

Ref: FORM-0630 Rev.16 9930042-E

ABX Minotrol 16



LOT MX 032 H
 Rev 1

CONTROL H

(Exp.) 2012-06-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.0	20.7							± 1.6
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	5.56	5.68							± 0.20
HB HGB	g/dl	17.8	18.0							± 0.6
	g/l	178	180							± 6
	mmol/l	11.05	11.18							± 0.37
HT HCT	%	50.6	50.6							± 3.0
	l/l	0.506	0.506							± 0.030
VGM MCV	µm ³ ; fl	91	89							± 4
TGMH MCH	pg	32.0	31.7							± 2.0
	f mol	1.99	1.97							± 0.12
CCMH MCHC	g/dl	35.2	35.6							± 3.0
	g/l	352	356							± 30
	mmol/l	21.85	22.08							± 1.86
IDR RDW	%	N/A	13.1							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	505	485							± 45
VPM MPV	µm ³ ; fl	N/A	7.2							± 2.0
LYMPHO	#	N/A	2.42							± 1.00
	%	N/A	11.7							± 5.0
MONO	#	N/A	1.08							± 0.80
	%	N/A	5.2							± 4.0
GRANULO	#	N/A	17.17							± 1.40
	%	N/A	83.0							± 7.0

Ref: FORM-0830 Rev.16 9930043-E