



ABX Minotrol 16



LOT **MX 111 L**
 Rev 1

CONTROL **L**

(Exp.) **2012-02-05**
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	1.8	1.9							± 0.4
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	2.44	2.49							± 0.15
HB HGB	g/dl	6.0	6.2							± 0.4
	g/l	60	62							± 4
	mmol/l	3.73	3.83							± 0.25
HT HCT	%	17.1	17.1							± 2.0
	l/l	0.171	0.171							± 0.020
VGM MCV	µm ³ ; fl	70	69							± 4
TGMH MCH	pg	24.6	24.7							± 2.0
	f mol	1.53	1.53							± 0.12
CCMH MCHC	g/dl	35.1	36.0							± 3.0
	g/l	351	360							± 30
	mmol/l	21.81	22.38							± 1.86
IDR RDW	%	N/A	14.6							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	72	88							± 20
VPM MPV	µm ³ ; fl	N/A	7.1							± 2.0
LYMPHO	#	N/A	0.97							± 0.40
	%	N/A	51.4							± 9.0
MONO	#	N/A	0.24							± 0.20
	%	N/A	13.0							± 7.0
GRANULO	#	N/A	0.66							± 0.40
	%	N/A	35.0							± 8.0

Ref: FORM-0830 Rev.16 9930041-E

ABX Minotrol 16



LOT MX 111 N
 Rev 1

CONTROL N

(Exp.) 2012-02-05
 (YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	7.3	7.5							± 0.8
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	4.68	4.78							± 0.18
HB HGB	g/dl	13.6	13.2							± 0.5
	g/l	136	132							± 5
	mmol/l	8.45	8.22							± 0.31
HT HCT	%	37.9	37.9							± 2.5
	l/l	0.379	0.379							± 0.025
VGM MCV	µm ³ ; fl	81	79							± 4
TGMH MCH	pg	29.1	27.7							± 2.0
	f mol	1.80	1.72							± 0.12
CCMH MCHC	g/dl	35.9	34.9							± 3.0
	g/l	359	349							± 30
	mmol/l	22.28	21.66							± 1.86
IDR RDW	%	N/A	13.8							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	250	253							± 40
VPM MPV	µm ³ ; fl	N/A	7.3							± 2.0
LYMPHO	#	N/A	1.99							± 0.60
	%	N/A	26.4							± 6.0
MONO	#	N/A	0.66							± 0.40
	%	N/A	8.8							± 5.0
GRANULO	#	N/A	4.88							± 0.60
	%	N/A	64.8							± 7.0

Ref: FORM-0830 Rev.16 9930042-E



ABX Minotrol 16



LOT **MX 111 H**
Rev 1

CONTROL **H**

(Exp.) **2012-02-05**
(YYYY - MM - DD)

		ABX VETPACK								
PARAMETRES PARAMETERS	UNITES UNITS	MICROS								TOLERANCES TOLERANCE
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ; 10 ⁹ /l	19.0	20.6							± 1.6
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	5.61								± 0.20
HB HGB	g/dl	18.0								± 0.6
	g/l	180								± 6
	mmol/l	11.18	11.18							± 0.37
HT HCT	%	50.5	50.5							± 3.0
	l/l	0.505	0.505							± 0.030
VGM MCV	µm ³ ; fl	90	88							± 4
TGMH MCH	pg	32.1	31.4							± 2.0
	fmol	1.99	1.95							± 0.12
CCMH MCHC	g/dl	35.7	35.6							± 3.0
	g/l	357	356							± 30
	mmol/l	22.14	22.13							± 1.86
IDR RDW	%	N/A	13.1							± 3.0
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	505	490							± 45
VPM MPV	µm ³ ; fl	N/A	7.6							± 2.0
LYMPHO	#	N/A	2.41							± 1.00
	%	N/A	11.7							± 5.0
MONO	#	N/A	0.86							± 0.80
	%	N/A	4.2							± 4.0
GRANULO	#	N/A	17.30							± 1.40
	%	N/A	84.0							± 7.0

Ref: FORM-0830 Rev.16 9930043-E